



Product Specification

DENARASE, 25 kU

Art. No.: 20804-25k

Version 1.3, valid as of October 05, 2020

Recombinant *Serratia marcescens* endonuclease produced by microbial fermentation with *Bacillus* sp. The production strain employed in the manufacturing of the product is a Genetically Modified Organism (GMO) of safety level S1.

The enzyme is supplied as liquid and formulated in 20 mM Tris-HCl pH 8.0 ±0.2, 20 mM NaCl, 2 mM MgCl₂, 50 % glycerol (v/v).

The product is produced under GMP conditions acc. to EU GMP regulations and filled under ISO 9001 conditions.

Parameter	Method	Specification
Appearance	visual	Clear, transparent solution
Activity	photometric ¹	> 250 U/μl
Purity	Protein purity determined by SDS-PAGE and silver staining	≥ 99 %
Protease activity	Protease detection assay	No protease activity detectable
Endotoxin level	LAL-Test acc. to Ph. Eur. 2.6.14, Method C	< 0.25 EU/kU
Total microbial count	TAMC/TYMC acc. to Ph. Eur. 2.6.12	Aerobic bacteria: < 5 cfu/200 μl Yeast/moulds: < 5 cfu/200 μl

Storage: Store at -20 °C.

Stability: Stable within specification range for a period of at least 24 months from the date of product release under proper storage conditions.

BSE / TSE / Animal derived material: The manufacturing process is free of materials with TSE/BSE risk and of raw materials from animal origin.

GMO-Statement: The product is free of the production strain.

Antibiotics: No antibiotics are used in the manufacturing process.

¹ Unit-Definition: One unit (U) will digest salmon sperm DNA to acid-soluble oligonucleotides equivalent to a ΔA260nm of 1.0 in 30 min at pH 8.0 at 37 °C.